

EEC Conformity Declaration

in accordance with
Annex III of the EEC Directive 73/23/EEC
and
Annex I of the EEC Directive 89/336/EEC

The company

Gerätebau Felix Schulte GmbH & Co. KG.
Rangestr. 46-48
D - 59581 Warstein

explained, that the products

Designation: Not motorized viewers:

EC#, EH#,XU#, XH#, XUH#, MOD#, MOD-Combined#, BW I#, BW II#, MAMMO#,
JX#46x48,
Iris#100 with switch, Iris#70 with switch, Iris#100 CWS with switch,

Motorized viewers:

ROTOLUX 120, ROTOLUX 120-8, ROTOLUX 220, ROTOLUX 220-8, ROTOLUX 320,
ROTOLUX 320-8, ROTOLUX 400, ROTOLUX 400-8,

correspond to the regulations of the CE guidelines specified above.

The products are in dependence of the harmonized standards:

- EN 60598-1:2005-03 Luminaires — Part 1: General requirements and tests
- EN 60598-2-25:2005-07 Luminaires — Part 2-25: Particular requirements — Luminaires for use in clinical areas of hospitals and health care buildings
- DIN EN 55015:2003-09 Limits and methods of measurement of radio disturbance characteristics of electrical lighting and similar equipment
- DIN EN 61000-3-2:2006-10 Electromagnetic compatibility (EMC) - Part 3-2: Limits - Limits for harmonic current emissions (equipment input current <= 16 A per phase)
- DIN EN 61000-3-3:2006-06 Electromagnetic compatibility (EMC) - Part 3-3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current <=16 A per phase and not subjected to conditional connection
- DIN EN 61547:2001-06 Equipment for general lighting purposes - EMC immunity requirements

Classification according to EN 61000-3-2:2006-06

All products except Iris#70 & Iris#100: Class C
Iris#70 & Iris#100: Class A

Gerätebau Felix Schulte
GmbH & Co. KG



ppa

Warstein, July 23, 2007