





PLANILUX®-Magnifiers for the reading of small details on radiographs

	Description	Order No																																							
	Hand magnifier, 120 mm diameter, twofold magnification with supersharp biconvex lens of optical glass with black contoured plastic handle	401004																																							
	Hand magnifier APLANAT, 65 mm diameter, fourfold with 2 lenses of optical glass with black contoured plastic handle especially suitable for use in mammography	401020																																							
	Hand magnifier, 55 mm diameter, sixfold <u>Light lens,</u> with black contoured plastic handle	401012																																							
	Small chain approx. 1,5 m long, with holder for magnifiers as specified above	401080																																							
	Magnifier mounted on a weight-compensated swivel arm designed for securing to the front of the film viewer or for wall-mounting, for the rest as specified under order No 401004 <u>This magnifier covers a viewing area width of 80 cm as a maximum!</u>	402008																																							
	Sliding focusable magnifier runs in a <u>ball-bearing roller guide</u> , can be moved across the <u>entire viewing area width</u> , <u>adjustable in height</u> , so that any position within the viewing area can be reached, lens as specified under order No 401004, available for the following film viewers:	<table border="1" data-bbox="564 1288 1329 1713"> <tbody> <tr> <td>Series D., E., X.. and combined viewers</td> <td>40x43 cm</td> <td>403401</td> </tr> <tr> <td></td> <td>80x43 cm</td> <td>403402</td> </tr> <tr> <td></td> <td>105x43 cm</td> <td>403403</td> </tr> <tr> <td></td> <td>120x43 cm</td> <td>403404</td> </tr> <tr> <td></td> <td>140x43 cm</td> <td>403405</td> </tr> <tr> <td></td> <td>43x105 cm</td> <td>403406</td> </tr> <tr> <td>Series DXHM</td> <td>72x46 cm</td> <td>403407</td> </tr> <tr> <td></td> <td>105x46 cm</td> <td>403411</td> </tr> <tr> <td></td> <td>140x58 cm</td> <td>403415</td> </tr> <tr> <td></td> <td>48x60 cm</td> <td>403409</td> </tr> <tr> <td></td> <td>108x60 cm</td> <td>403412</td> </tr> <tr> <td></td> <td>140x60 cm</td> <td>403416</td> </tr> <tr> <td>For other viewers:</td> <td colspan="2">Upon request!</td> </tr> </tbody> </table>	Series D., E., X.. and combined viewers	40x43 cm	403401		80x43 cm	403402		105x43 cm	403403		120x43 cm	403404		140x43 cm	403405		43x105 cm	403406	Series DXHM	72x46 cm	403407		105x46 cm	403411		140x58 cm	403415		48x60 cm	403409		108x60 cm	403412		140x60 cm	403416	For other viewers:	Upon request!	
Series D., E., X.. and combined viewers	40x43 cm	403401																																							
	80x43 cm	403402																																							
	105x43 cm	403403																																							
	120x43 cm	403404																																							
	140x43 cm	403405																																							
	43x105 cm	403406																																							
Series DXHM	72x46 cm	403407																																							
	105x46 cm	403411																																							
	140x58 cm	403415																																							
	48x60 cm	403409																																							
	108x60 cm	403412																																							
	140x60 cm	403416																																							
For other viewers:	Upon request!																																								
	Antiglare ring Of aluminum, matt black finish. For all magnifiers with 120 mm diameter.	405004																																							