



PLANILUX® - Mobile trolleys

	<p>Description</p> <p>Mobile trolley for film reading in <u>seated position</u> for the <u>rigid</u> mounting of a PLANILUX® viewer</p> <p>Dimensions in () for combined viewer with height of 1250 mm</p> <p>bolted square steel bar construction, with 4 castors, the two front castors are fitted with a brake, with storage plate and mounting instructions, for the following series/models:</p>
---	---

Series D., E., X..	Trolley external dimensions in mm			Top edge desk top in mm from floor top edge	Center viewing area height in mm (app.) from top floor edge		Order No
	W	H	D				
40x43 cm	600	888	715	720	1200		501201
80x43 cm	1000	888	715	720	1200		501202
105x43 cm	1250	888	715	720	1200		501203
120x43 cm	1400	888	715	720	1200		501204
140x43 cm	1600	888	715	720	1200		501205
Combined models					Lower unit	Upper unit	
2x 105x43 cm	1250	888	795	678 (628)	1060 (1013)	1662 (1638)	501206
2x 140x43 cm	1600	888	795	678 (628)	1060 (1013)	1662 (1638)	501208
DXHM							
72x46 cm	970	850	715	720	1200		501209
105x46 cm	1250	850	715	720	1200		501203
140x58 cm	1600	850	715	720	1200		501213
48x60 cm	680	850	715	720	1200		201211
108x60 cm	1280	850	715	720	1200		501212
140x60 cm	1600	850	715	720	1200		501213

	<p>Mobile trolley for film reading in <u>seated position</u></p> <p>With <u>pivot</u> for optimum setting of the viewing angle of a PLANILUX® viewer</p> <p>For the rest identical to model for rigid viewer mount.</p> <p>Please order always viewer together with trolley to enable us to make the necessary adaptations at factory.</p> <p>Not available for combined viewers.</p> <p>For the following series/models:</p>
---	--

Series D., E., X..	Trolley external dimensions in mm			Top edge desk top in mm from floor top edge	Center viewing area height in mm (app.) from top floor edge		Order No
	W	H	D				
40x43 cm	780	1203	715	720	1200		502201
80x43 cm	1180	1203	715	720	1200		502202
105x43 cm	1430	1203	715	720	1200		502203
120x43 cm	1580	1203	715	720	1200		502204
140x43 cm	1780	1203	715	720	1200		502205
DXHM							
72x46 cm	1150	1200	715	720	1200		502206
105x46 cm	1480	1200	795	720	1200		502203
140x58 cm	1780	1200	795	720	1200		502210
48x60 cm	860	1200	795	720	1200		502208
108x60 cm	1460	1200	795	720	1200		502209
140x60 cm	1780	1200	795	720	1200		502210